

## SINGLE USE ACCESSORIES FOR: **THE LEG & ARM LIFT**




**LIFTING STRAP  
2006-50**

Supports: ankle, calf,  
knee.



**LIFTING STRAP  
2006-60**

Supports: arm, fingers,  
forefoot and ankle.



**LIFTING STRAP  
2006-70-1**

Wide strap for various  
positions.



**ADHESIVE VELCRO  
COIN 2006-65**

For forefoot & finger  
lifting. Non allergenic.

To view our **technical accessories**  
please visit our website.

**Do you need more information  
about the Leg & Arm Lift?**


**Please get in touch**



mail@opitek.dk | +45 39 66 16 44

www.opitek-international.dk

CE MDR Class 1



INTRODUCING:

## **THE LEG & ARM LIFT**

Prevent injuries when lifting patients' arms and  
legs in the operating theatre



# THE LEG AND ARM LIFT



The preparation of a surgical site on either the arms or legs of a patient often requires lifting, holding and lowering procedures. This job is normally undertaken by theatre staff, and can be physically demanding.

This unique electrical lift uses a hand control to lift the patient's leg or arm into the correct position. The limb can also be lowered using the hand control to ensure a smooth and safe process.

- The lift operates by battery, and can easily be re-charged over-night.
- The leg lift can accommodate a patient in a supine, prone or sitting position.
- Bilateral lifting of 2 legs is possible with the standard equipped "2 leg lifting arm".
- A height extender (32cm) is included to help lifting arms when the patient is in a sitting position (ex.: Beach Chair).



## SINGLE USE LIFTING STRAPS

Can be used to support ankle, calf, knee, arm, fingers & fore-foot etc. Skin friendly.



## THE LIFT

Included with the lift are: standard lifting arm, bilateral lifting arm & height extender (32 cm).



Do you want to see the Leg & Arm Lift in 3D?



Scan the QR-code with your phone and watch.